




**RED ONE**  
MEDICAL

ECAT #: SPE2DE-23-D-A008



# LifeFlow<sup>®</sup> PLUS

## Blood & Fluid Infuser



The LifeFlow<sup>®</sup> PLUS Blood & Fluid Infuser is a compact, handheld infuser used for rapid and controlled blood or fluid infusion. Lightweight and easy-to-use, the LifeFlow Plus is designed to rapidly stabilize a critical patient consistent with TCCC and DCR protocols. LifeFlow PLUS can deliver 1 unit of blood or 500ml of fluid in 2 minutes at the point of injury and during transport without batteries or sensitive electronics.

### LifeFlow Plus Blood & Fluid Infuser

P/N: FG5005 (CS/5)

NSN: 6515-01-703-0004



- Improves the ability to deliver effective TCCC & DCR
- Lightweight, compact and portable; travels with patient
- Ideal for use in austere environments and transport
- Infuses blood products up to 10x faster than traditional methods
- An experienced user can set up in approximately 1 minute
- Operates with one hand, freeing second hand for other tasks
- Compatible with select blood & fluid warmers
- Engineered, designed and assembled in the USA

Available through Red One Medical on:

DAPA #: SP020018H0066 | GSA #: 47QSWA19D0073 | ECAT #: SPE2DE-23-D-A008

P: 912.662.1443

customerservice@redonemedical.com

www.redonemedical.com