


# LifeFlow PLUS® Rapid Infuser



## Product Description

LifeFlow PLUS is a portable, easy-to-use device for rapid delivery of blood or crystalloid to critically ill patients. The device is designed to provide careful titration of volume to meet individualized patient needs, enabling rapid resuscitation followed by immediate reassessment to inform subsequent clinical decisions. LifeFlow PLUS allows delivery of 1 unit of blood or 500mL of fluid in less than 2 minutes, and can be used at the point of injury, and during transport, and within a medical facility.

## Details

- Part Number: FG5005
- ECAT: SPE2DE-23-D-A008
- NSN: 6515-01-703-0004
- Unit of Issue U/I: Case (5 each)
- Distributor: Red One Medical
- DAPA: SP020018H0066
- Cage Code: 71UG3
- Manufacturer: 410 Medical
- Shelf Life: 3 years
- Weight: 0.75 lb (each)
- Cube: 11" H x 8" W x 3" D (each)
- Single patient use
- Engineered, designed, and assembled in the USA

## Benefits & Attributes

- Improves the ability to deliver effective TCCC & DCR
- Infuses blood products up to 10x faster than traditional methods
- Lightweight, compact and portable; travels with patient
- Ideal for use in austere environments and transport
- An experienced user can set up in approximately 1 minute
- Operates with one hand, freeing second hand for other tasks
- Works effectively with QinFlow Warrior portable warmer

# LifeFlow PLUS Rapid Infuser Replacement Tubing



## Product Description

The LifeFlow PLUS Blood & Fluid Replacement Tubing set is usable only with the LifeFlow PLUS Infuser. After four or more units of blood the flow rate with the original tubing may start to slow. Using the replacement tubing set allows providers to continue to use the same LifeFlow PLUS infuser handle but swap in a new tubing set for higher flow rates. Note: the LifeFlow® PLUS Blood & Fluid Infuser is an easy-to-use portable device for treating critically ill patients in need of rapid and controlled infusion of blood or crystalloids.

## Details

- Part Number: FG5006
- NSN: 6515-01-703-0024
- Unit of Issue U/I: Case (5 each)
- Distributor: Red One Medical
- Cage Code: 71UG3 ECAT, DAPA, PV (pending)
- Manufacturer: 410 Medical
- Shelf Life: 3 years
- Weight: 8 oz (each)
- Cube: 11" H x 8" W x 1.5" D (each)
- Single patient use
- Engineered, designed, and assembled in the USA

## Benefits & Attributes

- Improve TCCC (Tactical Combat Casualty Care) & DCR (Damage Control Resuscitation)
- Providers can infuse a unit of blood in 2 minutes, 3-10x faster than traditional methods
- Lightweight, compact and portable; travels with patient
- Ideal for use in austere environments and transport
- Easy-to-set up and use; a new user can set up in < 3 minutes
- Operate with one hand, freeing second hand for other tasks